


EXPLORE
FR α pha
Challenging FR α Staining Patterns

For more information, visit:
[Explore-FRalpha.com](https://explore-FRalpha.com)



This resource has been developed by AbbVie and is for illustrative purposes only. Images are representative and not intended for clinical use. Images are provided by Roche Diagnostics.




Heterogeneous FR α Expression

Characterization

Variable membrane staining intensity

Assessment Impact

May require closer examination at both higher and lower power to determine TC percentages




Obscuring Cytoplasmic Staining

Characterization

Staining in the cytoplasm of TCs

Assessment Impact

Challenging to determine membrane specific staining. Viewing the slide at a higher magnification may help elucidate staining patterns



Dot-like Staining Pattern


Characterization

Cases exhibiting dark brown secretions filling the gland lumens, mimicking membrane staining

Assessment Impact

Assessment for dot-like staining may be required for accurate interpretation. If present, dot-like staining of gland lumens is included in the FR α assessment. Dots observed **in the cytoplasm** are not included in assessment

This resource has been developed by AbbVie and is for illustrative purposes only. Images are representative and not intended for clinical use.




Apical Staining Pattern

Characterization

Staining in tissues with glandular or tubular morphology

Assessment Impact

If present, assessment of concurrent membranous staining in the portion of cells oriented towards the lumen of a gland or tubule is required



Non-specific Background Staining

Characterization

Staining not associated with FR α expression

Assessment Impact

Evaluation of internal control areas and positive/negative control slides should be used to assess the level of non-specific background staining

Tissue or Staining Artifacts

Characterization

Can include: fixation gradients and edge effects, DAB trapping, nuclear bubbling, lack of staining in regions, tearing or folding of tissue, loss of tissue section

Assessment Impact

Repeat staining or new sample may be required

DAB=3,3'-diaminobenzidine; FR α =folate receptor alpha.

Reference: VENTANA FOLR1 (FOLR1-2.1) Rx Dx Assay Interpretation Guide for Epithelial Ovarian, Fallopian Tube, or Primary Peritoneal Cancer. 1015089US Rev B. 2023-09-01.

© 2025 AbbVie Inc. All rights reserved.

ABBV-US-02008-MC, v1.0
Approved: April 2025